RacOOon Plate System

INDICAZIONI

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GB RacOOon Plate System

INDICAZIONI

The components of the RacOOon Plate System are used for surgical procedures of pediatric epiphysiodesis. The system drives epiphyseal growth, allowing for both angular correction (asymmetrical epiphysiodesis), and longitudinal correction (symmetrical epiphysiodesis), in case of deformities of the lower limbs, with a minimally invasive approach.

CONTROINDICATIONS

The RacOOon Plate System is a safe, standard, simple and reversible system. No diagnosis is a contraindication to its use, except for the ones reported below:

- Closed physis or less than 12 months of residual growth;
- Physiological axial deformity;
- Blount's disease, classified more than Langenskiold IV;
- Osteophytes causing the deformity (results of epiphysial detachment, etc.).

WARNINGS

The components of the RacOOon Plate System are single use, i.e. they are to be used only once. After explants, they must be eliminated. The surgeons and other users need to be aware that, if explanted, the retrieved devices must be treated as contaminated and contaminating products and that they need to be disposed according to the applicable procedures available at the hospital. The repeated use of any device leads to a system burdened by a previous use, which impairs the fatigue strength determined in the project and thus increases the risk of rupture. Gruppo Bioimpianti s.r.l. absolutely does not respond to any defects caused by such misuse.

IDENTIFICATION AND TRACEABILITY

All products manufactured by Gruppo Bioimpianti s.r.l. are identified by a label on the packaging which states the contents, size or similar, code, production batch, and sterility state. For devices supplied in a sterile state, the expiration date is also reported on the label. To guarantee external traceability even after implant surgery a series of labels are provided inside the packaging, to be used in the register of surgical operations and in patients' medical records.

PRESERVATION

Our devices should be kept in hygienically controlled surroundings in which the packaging will not be damaged. Should damage occur, the product should be returned, after due consideration, to Gruppo Bioimpianti S.r.l. for any necessary re-processing or elimination. The packaging system, however, guarantees the necessary protection against scratches, damages or loss of sterility, of the products supplied sterile, during routine handling. Do not use the device if the packaging is opened or damaged or in case the expiry date has passed, for sterile devices.

STERILISATION

The RacOOon Plate System components can be supplied by Gruppo Bioimpianti S.r.l. as sterile (according to the definition of the standard UNI EN 556) and as non sterile, depending on hospital needs. Different indications obviously need to be followed in the two cases.

Non-sterile product: the devices must be sterilised before use as per normal hospital practice (steam) and according to the normative UNI EN ISO 17665-1. If the product is not used during the surgery, it may be sterilised at a later date, provided it is handled with due care in the meantime and properly preserved (clean packaging in hygienically controlled surroundings). Should the product be impaired or show signs of superficial alteration, it should not be used.

Sterile product: products are supplied having undergone beta or gamma ray sterilisation according to the normative UNI EN ISO 11137-1,-2,-3. Sterility is indicated on the label together with the product's use-by date. The packaging (double pack) allows products to be handled by nonsterile personnel and the final stages of packaging elimination to be carried out by sterile operating personnel. If at any time during unpacking an anomalous situation arises, it is advisable not to use the piece because there is a risk that it may have lost the sterility. In this case the device has to be returned to Gruppo Bioimpianti that will carry out, if it is possible, a reconditioning.

MATERIALS

The components of the RacOOon Plate System are manufactured in biocompatible materials. The system is available in two versions: stainless steel (ISO 5832-1) and titanium alloy Ti6Al4V (ISO 5832-3). Our production process is subject to quality control appropriate to the sector (ISO 9001 e ISO 13485).

WARNINGS

To optimise results in the use of our products we draw attention to the following recommendations:

- handle products with due care in order to avoid accidental scratches, also during the operation itself;
- to fix the plate only use screws of the same material: do not combine screws in stainless steel with titanium plates and vice versa;
- always use two screws for each plate;
- it is responsibility of the surgeon to verify that the patient is not sensitive to the materials used to manufacture the implant's components;
- the components of the RacOOon Plate System should be removed when the pathology has been cured, subject to the patient's overall condition and in accordance with the physician decision;
- never combine components made by different manufacturers because there may be different tolerances that might prejudice the outcome of the implant;
- the devices manufactured by Gruppo Bioimpianti have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artefact in the MR environment. Their safety in the MR environment is unknown. Scanning a patient who has such devices may result in patient injury.

In general, it is always advisable to carry out a thorough pre-operative analysis to determine the most appropriate system and the right size compatible with the patient's clinical conditions.

Regarding surgical technique, Gruppo Bioimpianti S.r.l. provides indications on the use of its products which are different from the systems which have become part of traditional medical culture. The surgeon is to be thoroughly familiar with the implants, instruments and surgical procedures prior to performing surgery. In addition to product documentation, Gruppo Bioimpianti is willing to support the exchange of professional opinions between doctors who have already used such devices and those considering using them for the first time.

POST-OPERATIVE INDICATIONS

Full-load walking is allowed on the first post-operative day. Joint mobilization is immediate. If the patient feels insecure, walking with anti-brachial support sticks is recommended.

Post implant monitoring is recommended approximately every 3 months, although this activity remains at the discretion of the surgeon. The healing time for patients with idiopathic axial deviation is on average one calendar year, while for patients with dystrophic physis the healing times are difficult to quantify. Once the correct re-alignment of the mechanical axis of the limbs has been achieved it is necessary to remove the implant to avoid over correction.

In the event of a recurrence, the RacOOon Plate System can be used again by the same patient.

SURGICAL INSTRUMENTS

The instruments realized by Gruppo Bioimpianti are designed and manufactured to provide surgeons with the necessary instruments to realize a simple and safe implantation. Such instruments are supplied either as part of an operating set in dedicated trays or individually, and are identified by a code and a lot number, reported on the label when present and on the device if possible. Surgical instruments must be kept in hygiene-controlled premises, if possible in dedicated conservation units (surgical instrument trays). The materials used for manufacturing are biocompatible and guarantee good resistance to corrosion if stored in appropriate conditions. It should be remembered that organic fluids can cause local corrosion, so it is good practice to clean instruments after use, removing traces of blood or other organic fluids which might prove aggressive to the material if left in such conditions. The instruments, whether part of an operating set or supplied individually, are supplied non sterile. Gruppo Bioimpianti S.r.l. performs on each instrument a washing cycle according to the applicable Operative Instructions. Gruppo Bioimpianti underlines that it is responsibility of the hospital facility to perform on each instrument, before each surgery, the following sequence of processes, that must be performed appropriately according to the internal procedures approved by each facility: decontamination, washing (manual or ultrasound), re-washing, drying, control of absence of residual dirt or other impurities, control of lumen/holes, control of the integrity of the instruments, and sterilization according to a sterilization process conforming to the relevant standards (EN 556, ISO 11137-1,-2,-3 and ISO 11135-1 and UNI EN ISO 17665-1). Gruppo Bioimpianti S.r.l. manufactures instruments that can be sterilized several times, even using the simple systems generally found in hospitals (steam sterilizers). Should the instrument be impaired or show signs of superficial alteration, it should not be used. Any fracture or breaking of instruments shall be reported to the company. The use of components and/or instruments not included in the instruments set provided by Gruppo Bioimpianti S.r.l., except the generic instruments present in the operating rooms (drill, oscillating blades ...) could determine negative events. The instruments, if not carefully operated during the surgery, could cause damages to the implants. Surgical instruments are inevitably subject to wear with normal use and if they have experienced excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose, as illustrated in the surgical technique provided by Gruppo Bioimpianti S.r.l. It is recommended that all instruments be regularly inspected to identify any possible damage. **WARNING:** if the instruments are intended to be connected with an active medical device (electrical engine), the compatibility between two elements of the connection must be verified. Do not leave in the surgical site any trial component, packaging, instrument or part of it. Being medical devices, the instruments must be strictly used by medical personnel, which must be sufficiently familiar with their use before performing the surgery.