## MRS

### Modular Rotatory Shoulder Prosthesis

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
## IMPORTANT MEDICAL INFORMATION

## DEFINITIONS

Symbols may be used on the external label and/or on the package (UNI CEI EN ISO 15223-1). The following table provides the definition of these symbols.

Table 1: definitions of the possible symbols present on the label

Symbol	Definition
	Batch code
	Catalogue number
	Do not reuse
	Consult operating instructions
	Use by
	Do not reuse if package is damaged
	Date of manufacture
	Manufacturer
	Sterilized using ethylene oxide
	Sterilized using radiation

## DESCRIPTION

Modular Rotatory Shoulder endoprosthesis MRS provides to the surgeon a full and varied combinations with different dimensional characteristics for the treatment of tumor diseases affecting shoulder and / or its structures. It is composed by five elements: the humeral stem, the cylindrical spacer, the axial fixing screw, the bush and the spherical head.

## MATERIALS

Materials used for the manufacturing of Gruppo Bioimpianti products are conform to international standards of the industry, ensuring that, to date, the best levels of biocompatibility (ISO 5832). The production cycle of our equipment is complemented by a series of checks in accordance with industry standards (ISO 9001, ISO 13485). The materials used are in particular the following:

Table 2: materials used to manufacture the various elements of the MRS system.

Component	Material	Reference
Humeral stem	Ti6Al4V alloy	ISO 5832-3
Cylindrical spacer	Ti6Al4V alloy	ISO 5832-3
Axial fixing screw	Ti6Al4V alloy	ISO 5832-3
Spherical head	UHMWPE	ISO 5834-2
Bush	UHMWPE	ISO 5834-2

Note: all indicated standards are referred to the last current edition

## STERILIZATION

Product is supplied sterile (according to the definition reported on the standard UNI EN 556). In particular metal components are sterilized through rays (Beta or Gamma) according to the normative UNI EN ISO 11137-1;-2;-3 while the elements manufactured in UHMWPE are sterilized with ethylene oxide (ETO) according the normative UNI EN ISO 11135-1. Both sterilizations are guaranteed for 10 years. The sterilization method is indicated on the label in which is also reported the expiry date that indicates the

maximum time to use the product. The packaging (double pack) allows products to be handled by non-sterile personnel during conservation. The extraction of the product must be carried out by sterile operating personnel. If at any time during unpacking an anomalous situation arises, it is advisable not to use the piece because there is a risk that it may have lost the sterility and its return to Gruppo Bioimpianti which will carry out, if it is possible, a reconditioner.

#### **IDENTIFICATION AND TRACEABILITY**

All products manufactured by Gruppo Bioimpianti s.r.l. are identified by the label on the packaging which states the contents, size or another parameter which indicates the dimensions, the product's code, the production batch, the sterility symbol with the expiration date and the sterilization method. To guarantee external traceability even after implantation they are provided six labels inside the packaging, which shall be used on the surgery register and on the patients' medical records.

#### **PRESERVATION**

Our devices should be kept in hygienically controlled surroundings. The packaging system, however, guarantees the necessary protection against scratches, damages or loss of sterility of the product during routine handling. The handling and the preserving of the product have to be in order to not damage the packaging: If damage should occur, the product should be returned, after due consideration, to Gruppo Bioimpianti s.r.l. for any necessary re-processing or elimination. Do not use the product if the packaging is opened or damaged or in case the expiry date is passed.

#### **INDICATIONS**

The Modular Rotatory Shoulder Endoprosthesis MRS is used for conservative treatment of cancer. It allows shoulder replacement, without the application of the glenoid component, and a part of the humeral proximal diaphysis. The indications are in particular the following:

- Primary malignant tumors of the proximal humerus or scapula
- Aggressive benign tumors of the proximal humerus or scapula
- Secondary tumors with or without pathological fracture
- Malignant neoplasm of soft tissues in the region of the shoulder whose wide excision also requires the removal of a part or in total of the shoulder joint.

The humeral stem has to be fixed to the bone by cementation with bone cement (PMMA).

#### **CONTRAINDICATIONS**

The shoulder replacement is generally contraindicated in cases of:

- Acute or chronic infection, local or systemic,
- Serious deficiencies muscular, neurological or vascular of the joint
- Destruction or poor bone quality that may affect the stability of the system
- Allergies to bone cement, titanium alloy or polyethylene

#### **WARNINGS**

Never combine components made by different manufacturers because there may be different tolerances that might prejudice the outcome of the implant.

1. Handle the product with due care to avoid accidental scratches, even during surgery. An inadequate handling of the device, pre-operative and / or intra-operative, or any damage can lead to corrosion, fatigue failure or excessive wear. Use clean gloves when handling the device;
2. Use screws or accessories made of suitable material: in the case of a system composed of several elements, always use the same kit elements;
3. Follow the surgical technique provided by Gruppo Bioimpianti
4. A complete cleaning and removal of metal particles, bone cement, as well as any other debris resulting from surgery, it is necessary to minimize the wear of the articular surface;
5. Do not use the device if there is any damage;
6. Do not alter the system in any way;
7. Inadequate choice, positioning, alignment and fixation of one or more components could result in an incorrect distribution of the stresses, with a consequent reduction of life of the prosthetic implant.
8. An incorrect alignment or a plant not accurate can lead to a phenomenon of excessive wear and / or failure of the system itself.
9. All implantable devices of the MRS prosthesis are DISPOSABLE so it's not allowed their reuse. The repeated use leads to a system burdened by a previous use, which affects the fatigue strength and thus increases the risk of breakage;
10. Not implant a component in a patient if it had already been implanted, albeit temporarily, in another patient;
11. Use ALWAYS the appropriate locking screw for the MRS Prosthesis
12. Do not leave in the surgical site any trial, packaging, instrument or its part
13. It is the responsibility of the operating surgeon to determine whether there is adequate initial fixation and stability
14. Patient smoking may result in a delayed healing, non-healing and/or compromised stability;
15. It is the responsibility of the operating surgeon to ensure that patient is not sensitive to the materials used to manufacture the implant's components.
16. The surgeons and other users need to be aware that, if explanted, the retrieved devices must be treated as contaminated and contaminating products and that they need to be disposed according to the applicable procedures available at the hospital.
17. The devices manufactured by Gruppo Bioimpianti have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artefact in the MR environment. Their safety in the MR environment is unknown. Scanning a patient who has such devices may result in patient injury.

## **SURGICAL INSTRUMENTS**

The instruments sets realized by Gruppo Bioimpianti are designed and manufactured to provide surgeons with the necessary instruments to realize a simple and safe implantation. Such instruments are supplied either as part of an operating set in dedicated trays or individually, and are identified by a code and a lot number, reported on the label when present and on the device if possible. Surgical instruments must be kept in hygiene-controlled premises, if possible in dedicated conservation units (surgical instrument trays). The materials used for manufacturing are biocompatible and guarantee good resistance to corrosion if stored in appropriate conditions. It should be remembered that organic fluids can cause local corrosion, so it is good practice to clean instruments after use, removing traces of blood or other organic fluids which might prove aggressive to the material if left in such conditions. The instruments, whether part of an operating set or supplied individually, are supplied non sterile. Gruppo Bioimpianti S.r.l. performs on each instrument a washing cycle according to the applicable Operative Instructions. Gruppo Bioimpianti

underlines that it is responsibility of the hospital facility to perform on each instrument, before each surgery, the following sequence of processes, that must be performed appropriately according to the internal procedures approved by each facility: decontamination, washing (manual or ultrasound), re-washing, drying, control of absence of residual dirt or other impurities, control of lumen/holes, control of the integrity of the instruments, and sterilization according to a sterilization process conforming to the relevant standards (EN 556, ISO 11137-1,-2,-3 and ISO 11135-1 and UNI EN ISO 17665-1). Gruppo Bioimpianti S.r.l. manufactures instruments that can be sterilized several times, even using the simple systems generally found in hospitals (steam sterilizers). Should the instrument be impaired or show signs of superficial alteration, it should not be used. Any fracture or breaking of instruments shall be reported to the company. The use of components and/or instruments not included in the instruments set provided by Gruppo Bioimpianti S.r.l., except the generic instruments present in the operating rooms (drill, oscillating blades..) could determine negative events. The instruments, if not carefully operated during the surgery, could cause damages to the implants. Surgical instruments are inevitably subject to wear with normal use and if they have experienced excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose, as illustrated in the surgical technique provided by Gruppo Bioimpianti S.r.l. It is recommended that all instruments be regularly inspected to identify any possible damage.

**WARNING:** if the instruments are intended to be connected with an active medical device (electrical engine), the compatibility between two elements of the connection must be verified. Do not leave in the surgical site any trial component, packaging, instrument or part of it. Being medical devices, the instruments must be strictly used by medical personnel, which must be sufficiently familiar with their use before performing the surgery.

#### **PREOPERATIVE PRECAUTIONS**

In general, it is always advisable to carry out a thorough pre-operative analysis to determine the most appropriate system and the right size compatible with the patient's clinical conditions (age, physical activity/work, weight, mental illness, osteoporosis, vascular pathologies). The surgeon is to be thoroughly familiar with the implants, instruments and surgical procedures prior to performing surgery. Regarding surgical technique, Gruppo Bioimpianti s.r.l. provides indications on the use of its products. In addition to product's documentation, Gruppo Bioimpianti is willing to support the exchange of professional opinions between doctors who have already used such devices and those interesting on its using for the first time. The patient is to be made aware and warned of general surgical risks and possible adverse effects. The patient is to be advised of the limitations of the reconstruction and the need for protection of the implants from a full load until adequate fixation and healing have occurred. Any joint replacement system, including the implant/bone interface, cannot be expected to withstand activity levels and loads as would normal healthy bone and will not be as strong, reliable, or durable as a natural human joint. The patient should not have unrealistic functional expectations regarding to activities or occupations such as running, lifting weights and walking for long periods. Excessive, unusual and/or awkward movement and/or activity, trauma, excessive weight and/or obesity, could be implicated with premature failure of the implant. The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, has a finite expected service life and may need to be replaced at some time in the future. The patient should be advised that any noise or unusual sensation should be reported to surgeon as it may indicate implant malfunction.

#### **POSTOPERATIVE PRECAUTION**

The patient should be visited periodically after surgery and he has to be warned about the limitations of the postoperative process. The patient has to follow postoperative instructions, including those arising from follow-up visits, and he has to be aware that a correct rehabilitation is essential for the implant success. In

fact, in restoring all the functions of the joint, collaboration is needed to overcome the pathological state and, as required by the bio-mechanical specifications, the device must be integrated in such a way as to facilitate the normal processes of bone growth. To do this, it's important that the patient returns to a controlled normal activity, through exams and radiological analysis, in order to manage the loads that can be applied to the joint and to the medical device. Except initially, our devices do not create any particular problems for the patient in everyday life. Should any problems arise, please inform immediately the patient. Periodic postoperative X-ray testing is recommended.

**POSSIBLE ADVERSE EFFECTS**

1. Thromboembolic
2. Cerebral Stroke
3. Myocardial Infarction
4. Mobilisation prosthetic component (glenoid or humeral stem)
5. Dissection of the prosthetic components
6. Glenoid erosions
7. Rupture of the rotator cuff
8. Intra-operative bone fractures particularly in the presence of low bone density
9. Fractures of the prosthetic implant
10. Dislocation or subluxation
11. Vascular and nerve injuries
12. Infection early or late superficial or deep
13. Ectopic ossification
14. Instability
15. Joint stiffness

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GRUPPO BIOIMPIANTI S.r.l.