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K-MONO *system* Unicompartimental femoro-tibial knee prosthesis

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IMPORTANT MEDICAL INFORMATION**K-MONO System**

Unicompartmental femoro-tibial knee prosthesis

DEFINITIONS

Symbols may be used on the external label and/or on the package (UNI CEI EN ISO 15223-1). The following table provides the definition of these symbols.

Table 1: definitions of the possible symbols present on the label

Symbol	Definition
	Batch code
	Catalogue number
	Do not reuse
	Consult operating instructions
	Use by
	Do not reuse if package is damaged
	Date of manufacture
	Manufacturer
	Sterilized using ethylene oxide
	Sterilized using radiation

DESCRIPTION**K-MONO:**

The unicompartmental prosthesis K-MONO with fixed articular bearing is a knee prosthesis for primary indication to be used in only one femoro-tibial compartment (lateral or medial) of the knee. Such prosthesis is composed by different components, such as the femoral component, the tibial tray and the fixed articular insert. The femoral and tibial components are available in different designs:

- Femoral component:
 - 3CUT
 - RES
- Tibial component:
 - MetalBack
 - AllPoly

In addition, a tibial plug, preassembled to the cementless MetalBack tibia, and a screw to be used with such version of the tibial tray, are also available.

La femoral component 3CUT and the tibial component MetalBack in Titanium are available in both the Cemented and Cementless versions, the latter with a porous coating. The femoral component RES, the

tibial component MetalBack in CrCoMo and the tibial component AllPoly are only available in the cemented version.

The femoral components 3CUT e RES are also available in the anallergic version, which is characterized by a TiNbN coating and meant to be cemented. Such coating makes the components manufactured in CoCrMo anallergic, while the components manufactured in Ti6Al4V are already anallergic even without the coating.

The femoral component 3CUT is available also in the uncemented anallergic version, characterized by the presence of both a TiNbN coating and a porous coating.

MATERIALS

The materials used to manufacture Gruppo Bioimpianti's products are compliant with the International Standards which guarantee, nowadays, the best Biocompatibility levels (Table 2). The production process of our products includes a series of controls which are conforming to the applicable standards (ISO 9001, ISO 13485).

Table 2: materials used to manufacture the various elements of the K-MONO.

Component	Material	Reference
Femoral component (3CUT and RES)	CoCrMo alloy	ISO 5832-4
Tibial tray Metal Back CrCoMo	CoCrMo alloy	ISO 5832-4
Tibial Tray Metal Back Titanium	Ti6Al4V alloy	ISO 5832-3
Tibial tray AllPoly	UHMWPE XLPE XLPE + Vitamina E	ISO 5834-2 ISO 5834-2 ISO 5834-2 / ASTM F2695
Insert	UHMWPE XLPE XLPE + Vitamina E	ISO 5834-2 ISO 5834-2 ISO 5834-2 / ASTM F2695
Tibial plug for Cementless Tibial tray Metal back	Ti6Al4V alloy	ISO 5832-3
Screw for Cementless Tibial tray Metal back	Ti6Al4V alloy	ISO 5832-3

STERILIZATION

Product is supplied sterile (according to the definition reported on the standard UNI EN 556). In particular metal components are sterilized through rays (Beta or Gamma) according to the normative UNI EN ISO 11137-1;-2;-3 while the elements manufactured in polyethylene are sterilized with ethylene oxide (ETO) according the normative UNI EN ISO 11135-1. Both sterilizations are guaranteed for 10 years. The sterilization method is indicated on the label in which is also reported the expiry date that indicates the maximum time to use the product. The packaging (double pack) allows products to be handled by non-sterile personnel during conservation. The extraction of the product must be carried out by sterile operating personnel. If at any time during unpacking an anomalous situation arises, it is advisable not to use the piece because there is a risk that it may have lost the sterility and its return to Gruppo Bioimpianti which will carry out, if it is possible, a reconditioning.

IDENTIFICATION AND TRACEABILITY

All products manufactured by Gruppo Bioimpianti s.r.l. are identified by the label on the packaging which states the contents, size or another parameter which indicates the dimensions, the product's code, the production batch, the sterility symbol with the expiration date and the sterilization method. To guarantee external traceability even after implantation they are provided a series of labels inside the packaging, which shall be used on the surgery register and on the patients' medical records.

PRESERVATION

Our devices should be kept in hygienically controlled surroundings. The packaging system, however, guarantees the necessary protection against scratches, damages or loss of sterility of the product during routine handling. The handling and the preserving of the product have to be in order to not damage the packaging: If damage should occur, the product should be returned, after due consideration, to Gruppo Bioimpianti S.r.l. for any necessary re-processing or elimination. Do not use the product if the packaging is opened or damaged or in case the expiry date is passed.

INDICATIONS

K-MONO:

The unicompartmental knee prosthesis with fixed bearing K-MONO is indicated for primary surgery for the unicompartmental femoro-tibial replacement of the knee joint. It is indicated in case of pathologies that affect only one of the femoro-tibial compartment of the knee, such as: unicompartmental arthrosis, post-traumatic unicompartmental arthrosis, localized necrosis limited to the femoral condyle.

Such prosthetic design must be used when all the ligaments of the knee joint involved are anatomically and functionally intact.

The femoral and tibial components which do not have porous coating are EXCLUSIVELY indicated for a cemented application. The metallic components with porous coating, instead, are indicated for a cementless application. In case of use of a cementless tibial tray with porous coating, the use of a self-threading cancellous screw is mandatory.

The anallergic femoral component is characterized by a TiNbN coating that makes it suitable for use in patient with sensibility to metal ions and is indicated for a cemented application. The 3CUT femoral component is also available in the uncemented anallergic version, for application without bone cement in patients that show sensibility both to metal ions and bone cement.

The cemented anallergic tibial component is represented by the tibial tray MetalBack in Titanium and the tibial tray AllPoly.

CONTRAINDICATIONS

Primary knee arthroplasty surgeries are absolutely contraindicated in case of: local or systemic infection, sepsis or osteomyelitis. It is relatively contraindicated in the following cases:

1. Osteoporosis;
2. Injuries of the bone structures (i.e., giant tumor cell or any malignant tumor, bone cysts, etc..);
3. Uncooperative patient or a patient with neurologic disorders, who is unable to follow instructions;
4. Systemic and/or metabolic disorder which lead a progressive deterioration of the bone support;
5. Osteomalacia;
6. Distant foci of infections which may spread to the implant site;
7. General Weakened of the immune system;
8. Vascular insufficiency, muscular atrophy or neuromuscular disease.;
9. Incomplete or deficient soft tissue surrounding the knee;
10. Usual and excessive consumption of alcohol and / or drugs;
11. Obesity;
12. Skeletal immaturity.

In addition to the abovementioned generic contraindications, the following contraindications, specific for the unicompartmental femoro-tibial knee prosthesis, also apply:

1. Adiposity;
2. Severe dysfunction of the two cruciate ligaments;
3. Previous corrective osteotomy;
4. Deficit of extension of more than 10 °;
5. Advanced pathology in the other compartments;
6. Rheumatoid arthritis;
7. Retro-patellar arthrosis;
8. Insufficient bone substance and poor bone quality which could compromise a stable anchoring of the implant;
9. Varus/valgus deformity of more than 10°.

WARNINGS

1. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions and a subsequent reduction in the service life of the prosthetic implant;
2. Malalignment of the components or an inadequate implant can place inordinate forces on the components which may cause excessive wear and/or the failure of the implant;
3. Malalignment or soft tissue imbalance can place inordinate forces on the components, which may cause excessive wear. In this case revision surgery may be required to prevent component failure;
4. Handle products with the utmost care to avoid scratching or damaging the implant, also during surgery. Improper preoperative or intraoperative implant handling or damage (scratches, dents, etc.) can lead to crevice corrosion, fretting, fatigue fracture and/or excessive wear. Use clean gloves when handling implants;
5. A complete cleaning and removal of metallic particles, bone cement as well as any surgical debris is necessary to minimize articular surface wear;
6. Do not use the implant if there is some damage;
7. DO NOT modify implants;
8. The K-MONO components are DISPOSABLE so it's not allowed their reuse. The repeated use leads to a system burdened by a previous use, which impairs the fatigue strength and thus increases the risk of rupture;
9. Do not treat patients with components which have been placed, even if momentarily, in a different patient;
10. DO NOT USE a K-MONO prosthetic component with a component manufactured by another company because compatibility is not ensured;
11. ALWAYS use the self-threading cancellous screw in case of implant of a cementless tibial tray;
12. The tibial tray MetalBack and the related articular insert MUST be of the same size;
13. There is a complete compatibility between the sizes of the tibial tray component (MetalBack or AllPoly) and the femoral component;
14. Particular care should be taken to ensure the complete support of all cemented components to avoid an overstressing which could lead to a procedure failure;
15. It is the responsibility of the operating surgeon to determine whether there is adequate initial fixation and stability.
16. Patient smoking may result in a delayed healing, non-healing and/or compromised stability;
17. It is the responsibility of the operating surgeon to ensure that patient is not sensitive to the materials used to manufacture the implant's components.

18. The surgeons and other users need to be aware that, if explanted, the retrieved devices must be treated as contaminated and contaminating products and that they need to be disposed according to the applicable procedures available at the hospital.
19. The devices manufactured by Gruppo Bioimpanti have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artefact in the MR environment. Their safety in the MR environment is unknown. Scanning a patient who has such devices may result in patient injury.

The K-MONO system provides the surgeons with the necessary instruments to reduce pain and restore the knee joint function. While these products are generally successful in attaining these goals, they cannot be expected to withstand the activity levels and function of a normal healthy joint.

SURGICAL INSTRUMENTS

The instruments sets realized by Gruppo Bioimpanti are designed and manufactured to provide surgeons with the necessary instruments to realize a simple and safe implantation. Such instruments are supplied either as part of an operating set in dedicated trays or individually, and are identified by a code and a lot number, reported on the label when present and on the device if possible. Surgical instruments must be kept in hygiene-controlled premises, if possible in dedicated conservation units (surgical instrument trays). The materials used for manufacturing are biocompatible and guarantee good resistance to corrosion if stored in appropriate conditions. It should be remembered that organic fluids can cause local corrosion, so it is good practice to clean instruments after use, removing traces of blood or other organic fluids which might prove aggressive to the material if left in such conditions. The instruments, whether part of an operating set or supplied individually, are supplied non sterile. Gruppo Bioimpanti S.r.l. performs on each instrument a washing cycle according to the applicable Operative Instructions. Gruppo Bioimpanti underlines that it is responsibility of the hospital facility to perform on each instrument, before each surgery, the following sequence of processes, that must be performed appropriately according to the internal procedures approved by each facility: decontamination, washing (manual or ultrasound), re-washing, drying, control of absence of residual dirt or other impurities, control of lumen/holes, control of the integrity of the instruments, and sterilization according to a sterilization process conforming to the relevant standards (EN 556, ISO 11137-1, -2, -3 and ISO 11135-1 and UNI EN ISO 17665-1). Gruppo Bioimpanti S.r.l. manufactures instruments that can be sterilized several times, even using the simple systems generally found in hospitals (steam sterilizers). Should the instrument be impaired or show signs of superficial alteration, it should not be used. Any fracture or breaking of instruments shall be reported to the company. The use of components and/or instruments not included in the instruments set provided by Gruppo Bioimpanti S.r.l., except the generic instruments present in the operating rooms (drill, oscillating blades..) could determine negative events. The instruments, if not carefully operated during the surgery, could cause damages to the implants. Surgical instruments are inevitably subject to wear with normal use and if they have experienced excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose, as illustrated in the surgical technique provided by Gruppo Bioimpanti S.r.l. It is recommended that all instruments be regularly inspected to identify any possible damage.

WARNING: if the instruments are intended to be connected with an active medical device (electrical engine), the compatibility between two elements of the connection must be verified. Do not leave in the surgical site any trial component, packaging, instrument or part of it. Being medical devices, the instruments must be strictly used by medical personnel, which must be sufficiently familiar with their use before performing the surgery.

PREOPERATIVE PRECAUTIONS

In general, it is always advisable to carry out a thorough pre-operative analysis to determine the most appropriate system and the right size compatible with the patient's clinical conditions (age, physical activity/work, weight, mental illness, osteoporosis, vascular pathologies). The surgeon is to be thoroughly familiar with the implants, instruments and surgical procedures prior to performing surgery. Regarding surgical technique, Gruppo Bioimpianti S.r.l. provides indications on the use of its products. In addition to product's documentation, Gruppo Bioimpianti is willing to support the exchange of professional opinions between doctors who have already used the K-MONO system and those interesting on its using for the first time.

It is important to realize an appropriate patients selection, during which the following important factors should be considered:

1. The patient must have reached full skeletal maturity;
2. Need to obtain pain relief and improve function;
3. Ability and willingness of the patient to follow indications, including control of weight and adequate activity level;
4. Patient's weight;
5. Activity and occupation level;
6. State of health of the patient.

The patient is to be made aware and warned of general surgical risks and possible adverse effects. The patient is to be advised of the limitations of the reconstruction and the need for protection of the implants from a full load until adequate fixation and healing have occurred. The patient should not have unrealistic functional expectations regarding to activities or occupations such as running, lifting weights and walking for long periods. Excessive, unusual and/or awkward movement and/or activity, trauma, excessive weight and/or obesity, could be implicated with premature failure of the implant.

POSTOPERATIVE PRECAUTION

The patient should be visited periodically after surgery and he has to be warned about the limitations of the postoperative process. The patient has to follow postoperative instructions, including those arising from follow-up visits, and he has to be aware that a correct rehabilitation is essential for the implant success. In fact, in restoring all the functions of the joint, collaboration is needed to overcome the pathological state and, as required by the bio-mechanical specifications, the device must be integrated in such a way as to facilitate the normal processes of bone growth.

To do this, it's important that the patient returns to a controlled normal activity, through exams and radiological analysis, in order to manage the loads that can be applied to the joint and to the medical device.

Except initially, our devices do not create any particular problems for the patient in everyday life. Should any problems arise, please inform immediately the patient.

Periodic postoperative X-ray testing is recommended.

POSSIBLE ADVERSE EFFECTS

The possible adverse effects of the K-MONO system are similar to those occurring with any total knee replacement, both revision and primary surgeries, and include the following:

1. General Complications related to surgical interventions, medications, other instruments used, blood etc.
2. Particulate wear debris from metallic and polyethylene components may be present in adjacent tissue or fluid. It has been reported that wear debris may initiate a cellular response resulting in osteolysis which can lead to the implant loosening;
3. Early or late postoperative infection and allergic reaction;
4. Intraoperative bone fracture may occur, particularly in the presence of poor bone stock caused by osteoporosis, bone defects from previous surgery, bone resorption, or while inserting the device;
5. Loosening, migration or fracture of the implants can occur due to trauma, loss of fixation, malalignment, mal position, bone resorption, unusual and/or awkward movement and/or excessive activity;
6. Periarticular calcification or ossification, with or without impediment of joint mobility;
7. Inadequate range of motion due to improper selection or positioning of components;
8. Variation in the leg length;
9. Dislocation or subluxation due to: inadequate fixation, malalignment or mal position, unusual and or excessive movement, trauma, obesity, weight gain. Muscle and soft tissue laxity;
10. Fatigue fracture of a component can occur as a result of loss of fixation, strenuous activity, trauma, malalignment, or excessive weight;
11. Corrosion at the interfaces between components;
12. Wear and/or deformation of articulating surfaces;
13. Varus-valgus deformity;
14. Patellar tendon rupture and ligamentous laxity;
15. Intraoperative or postoperative bone fracture, and/or postoperative pain.
16. Hematomas and delayed wound healing;
17. Temporary or permanent neuropathies;
18. Venous thrombosis and pulmonary embolism;
19. Cardiovascular disorders;
20. Occasional and permanent nerve dysfunction.

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GRUPPO BIOIMPIANTI S.r.l.