Enjoy Mobility

# EVO-NAIL PLUS

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## EVO-NAIL PLUS: Equal in appearance, New in operating room

EN

EVO-NAIL PLUS locked nail for proximal fractures (Standard, Intermediate or Long version) is an alternative to open synthesis techniques. It provides reduction and immobilization of the fracture by a closed technique which limits the risk of infection and blood loss. As is normal with all orthopaedic implant systems, there is a risk of failure and this can vary from patient to patient.

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PRE-ASSEMBLED LOCKING MECHANISM
SIMPLIFIED SURGICAL TECHNIQUE
RADIOLUCENT INSERTION HANDLE WITH RADIOGRAPHIC INDICATORS
LOW VALGUS ANGLE (4°)
STATIC OR DYNAMIC PROCEDURE
THREE IMPLANTS LENGTHS AVAILABLE

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## Indications

EN

Choosing this type of implant and the size to use is a matter of general traumatological principles and of the surgeon's experience. Generally speaking, the EVO-NAIL PLUS is used for:

- 1) BASICERVICAL FRACTURES
- 2) PERTROCHANERIC FRACTURES
- 3) INTERTROCHANERIC FRACTURES
- 4) SUB-TROCHANERIC FRACTURES

If any of these is accompanied by a diaphyseal fracture, the Evo-Nail Plus Long is used.

The implant should not be used when there is:

- inadequate quantity and quality of bone tissue
- infection in progress
- excessive geometrical anomaly in the shape of the femur

## Materials

EN

Evo-Nail Plus and all its components is made of titanium alloy Ti 6Al 4V grade 5 ELI according to ISO 5832/3. Surface treatment with anodization type II DOTIZE®. The pre-assembled locking mechanism, is equipped with a thin polyethylene UHMWPE band, which prevents the possible loosening.

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## Pre-operative planning and angle selection

EN

First of all choose the angle between the nail and the transcervical component (screw).

On the basis of X-rays of the fractured femur and the contra-lateral intact one and with the help of the pre-operative slides, the choice of angle between the endomedullary diaphyseal component (nail) and the transcervico-cephalic component must be either 125° or 130°. (FIG.1)

This choice must be made very carefully because it affects the final position of the synthesis system in the femur head and therefore of the strength of its hold and stability over time.

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## Access



EN

Introduce the nail, from the apex of the great trochanter through the common supratrochanteric access and then, under X-ray control, prepare an entry for the diaphyseal endo-midullary component (nail). (FIG.1)

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### Canal perforation

EN

Having decided on the point of introduction of the nail, use the Curved Cannulated Awl (140030641) (FIG.3) to perforate, checking both views with the image amplifier as you go.

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## Conical reamer sleeve positioning

EN

Remove the Awl and position the Conical Reamer Sleeve (140040470) (FIG.5) for the conical reamer.

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## Canal reaming


EN

Insert the Conical Reamer (140040467) on the guide wire still in place. Ream the endomedullary canal as far as the mark on the stem of the reamer. (FIG.6)

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## Nail assembly to the nail inserter

EN

Assemble the nail on the appropriate Nail Inserter (140040404 or 140040405) using the nail Fixing Bolt (140040413) (FIG.7) and the Ball Tip Screwdriver (140040425).

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## Nail insertion

EN

Manually introduce the nail fixed to its inserter into medullary cavity of the femur. Drive in the nail till the axis of the lag screw is positioned in the lower quadrant of the femur neck. (FIG.8)  
The nail should advance without difficulty by simply pushing with your hand.

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## Insertion with the help of the impactor (optional)

EN

Should this insertion require a slight force, use the Impactor (140040495), which must be placed on top of the nail fixing bolt. (FIG.9)

- Carry out and X-ray check.
- Avoid any kind of percussion against the Nail Inserter.

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### Lag screw guide sleeve insertion

EN

Remove the wire and Insert the lag Screw Guide Sleeve (140040430) in the Nail Inserter till it touches the bone. The lag Screw Guide Sleeve must stay in place fixed to the nail inserter with its threaded knob, until just before the introduction of the distal screws. It is important to introduce the Proximal Awl (140040440) and then penetrate the lateral cortical. (FIG.10)

- Fix the sleeve to the nail inserter using its threaded knob to stabilize the system.
- Don't tighten the threaded knob excessively as this may deform the sleeve.

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

### Optional Steinmann nail insertion

EN

To stop the femur head twisting when the lag screw is inserted, a 4.0 mm diameter Steinmann Nail (140040455 - L.285 mm) or (140040456 L.300 mm) may be put into the Nail Inserter using the Steinmann Guide Sleeve (140040450). (FIG.11)

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### Insertion of the cannulated proximal probe and threaded guide wire

EN

Remove the Awl and insert the Cannulated Proximal Probe (140040435) into the Guide Sleeve with the 3x410 mm Threaded Guide Wire (140040445). (FIG.12)

The probe should be pushed in till it touches the bone. Push in the wire through the lateral cortical till it reaches the subchondral zone of the femur head.

- It is recommended not press too much on the guide wire to prevent it to take wrong directions.
- If the X-rays show the guide wire is bent it must absolutely be removed and re-inserted correctly or even replaced with another one.
- Carry out X-ray checks (anteroposterior and axial).

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### Selection of screw lag length



EN

Having correctly positioned the guide wire, measure the part of it remaining outside the probe by placing it in the slot of the lag Screw Gauge (140040575) resting against the end of the probe itself. (FIG.13)

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### Lag screw housing creation

EN

Extract the probe from the Guide Sleeve being careful that the guide wire too doesn't move back. Introduce the Lag Screw Step Drill (140040475) with a ring on the back to vary the length of the perforating part of the reamer on the basis of the measurement taken (step 13). (FIG.14)

- If the bone is osteoporotic, it's possible to use the lag Screw Solid Drill (140040477).

Drill until coming into contact with the sleeve, thus obtaining a housing for the lag screw.

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### Lag screw insertion

EN

Having extracted the Lag Screw Drill, keeping the guide wire well in place, put the screw in through the Lag Screw Guide Sleeve. To do this, use the Lag Screwdriver (140040485). Screw it till the symbol marked on the stem (◆◆◆) of the screwdriver reaches the edge of the Sleeve. (FIG.15)  
Twist the ring Clockwise as far as the stop and proceed with compression when necessary.

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
### Locking of the lag screw

EN


Having inserted the lag screw, block it with the pre-assembled headless screw, using the Screwdriver (140040431).

- Make sure the lag screw is properly locked by trying to screw it. (FIG.16)

### Distal hole choice

EN

Insert the Distal Guide Sleeve (140030585) in the chosen hole in the nail inserter till it touches the cortical. "Round" for the round hole, "Stat" (static) for the proximal part of the slot and "Din" (dynamic) for distal part of the slot. If the 250 mm nail is used, choose the holes provided for that nail. (FIG.17)  
Introduce the Distal Awl (140030610) and make a hole in the cortical to invit the tip of the drill.

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### Creation of the hole for distal screw

EN

Remove the awl and put into the sleeve the Distal Guide Sleeve (140040460) with the Graduated Drill Bit (140040480) inside. (FIG.18)

- Perforate the two corticals advancing slowly.
- Carry out X-ray checks.

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### Measuring the depth of the hole


EN

Measure the depth of the hole for the distal screw housing using the graduated scale on the Drill Bit (140040479). (FIG.19)  
The reading of the measure has to be done with the tip positioned approximately 1 mm outside of the second cortical.  
The reading of the measure determines the length of the screw to use.

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## Distal screw insertion

EN

Having determined the distal screw length, introduce the screw using the Distal Screw Screwdriver (140040481) through the distal sleeve till the reference notch. (FIG.20)

## Optional insertion of the screw cap

EN

When the Nail Inserter is removed, the Nail Screw Cap (140400011) can be applied to the proximal end of the nail to protect the thread for eventual extraction. (FIG.21)

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

### EVO NAIL PLUS extraction

EN

If the EVO NAIL PLUS nail needs to be removed, unscrew the distal screws with the Distal "T" screw screwdriver (140040481).

Using the screwdriver (140040431), unscrew the cap (140400011) from the proximal end of the nail and unblock the Lag screw by unscrewing the pre-assembled headless screw.

Extract the lag screw using the Lag screwdriver (140040485) in extraction.

Once the Lag screw has been removed, it is possible to remove the nail: screw the Extraction bolt (140040415) to the proximal end and screw the Sliding hummer (140030665) onto the bolt.

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# EVO-NAIL PLUS

Vassoi porta strumenti  
Trays for instruments

IT  
EN  
ES

Bandeja  
para instrumentos


Vassoi porta strumenti  
Trays for instruments

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



Instrument set  
EVO NAIL PLUS EXTRACTION KIT  
REF: 140040523; 140040524

IT	EN	ES	REF
29 Massa battente	Sliding hummer	Martillo	140030665
30 Cacciavite a T Es.3.5	T Screwdriver Hex 3.5	Destornillador en T Hex 3.5	140040481
31 Cacciavite vite distale Es.5	T Screwdriver Hex 5	Destornillador en T Hex 5	140030615
32 Cacciavite grano di bloccaggio	Screwdriver pre-assembled screw	Destornillador tornillo pre-ensamblado	140040420
33 Inseritore vite cefalica	Lag screwdriver	Introductor	140040485
34 Chiave cardanica	Cardan key	Llave cardan	140030600
35 Bullone per estrazione	Extraction bolt	Bulón de extracción	140040415

Componente diafisaria standard (195 mm)	IT	
Standard Diaphyseal component (195 mm)	EN	
Componente diafisaria estándar (195 mm)	ES	
REF	CCD	DISTAL Ø
140810125	125°	10 mm
140810130	130°	10 mm
140810135*	135°	10 mm
140811125	125°	11 mm
140811130	130°	11 mm
140811135*	135°	11 mm

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## IFU

Use the QR-Code to view complete product information, including instructions for use, indications and contraindications, precautions and warnings

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This surgical technique is exclusively intended for medical professionals, especially physicians and surgeons.

This document does not constitute medical advice, it does not dispense medical recommendations and it does not convey any diagnostic or therapeutic information.

Informations and techniques presented in this document were compiled by a team of medical experts and Gruppo Bioimpianti's specialists; however Gruppo Bioimpianti excludes any liability for improper use of informations.

For any information or enquires about this publication or anything else, contact GRUPPO BIOIMPIANTI.

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