Enjoy Mobility

K-MOD™ **REV** KNEE-MODULAR REVISION

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The modular K-MOD REV system is designed for primary and revision knee replacement where there is bone loss or defects, correctable misalignment or instability, aseptic mobilization, osteolysis or infection. Thanks to its full modularity with K-MOD primary implants, the system can be used in a wide range of situations requiring lower invasiveness, bone saving, correct alignment, balancing of extension and flexion gaps and restoration of the knee's biomechanics.

K-MOD REV system allows to intraoperatively choose the most appropriate degree of constraint for the individual patient:

- Anterior and posterior stabilization using Ultra-Congruent UC inserts or Dynamic Congruence DC medial pivot inserts (BOXLess femur) or PS posterior stability inserts (CCK femur). In this version the system is intended for patients with collateral ligaments anatomically and functionally present.
- Constriction in varus/valgus and in femorotibial rotation using the CCK inserts combined with CCK femur. In this version, the system is intended for patients with collateral ligaments present but who need mechanical stabilization through the prosthetic implant.

Both versions use the same stems, offset modules, wedges, the same K-MOD REV tibial component and the same K-MOD patellar component.

The operating technique requests the same K-MOD REV femoral resections (anterior, posterior, distal, chamfers): proceeding sequentially, the surgeon will be able to evaluate the best degree of constraint, even in the final phase with the trial reduction. The transition to the CCK constraint requires only the preparation of the central box, that is performed using a special guide centered at the intra-medullary level.

After determining the degree of constraint, it is possible to perform the final verification of the most suitable insert, even with the definitive femoral and tibial components already implanted.

In case of K-MOD REV CCK implant, in order to optimize the final femorotibial rotation, an alternative procedure to the classic sequence used for the K-MOD REV is reported.

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On the patient's most recent X-rays perform thorough pre-op planning on the AP and ML planes by using the X-Ray templates: K-MOD REV BoxLess X-Ray Template Fem AP/ML, Tib AP/ML, Stem (KREV-XR-000/004).


Evaluate (magnification: 15%): sizes of the femoral and tibial components, diameters/lengths of the diaphyseal stems and whether off-set modules and femoral and/or tibial augmentations are needed.

The X-Ray templates have a useful highlighted area for the connection between components and stems or stems with a 3 or 6mm off-set module.

Evaluation is indicative and must always be checked and modified during surgery if necessary.

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To simplify the nurse's intra-op tasks and considering the full compatibility with the K-MOD primary system, the trays come in the following logical sequence.

Tibial implant only. Only the following trays are required:

- REAMERS (300113900)
- TIBIAL INSTRUMENTS (300113902)
- COMMON INSTRUMENTS (300113903)
- TRIAL STEMS (300113901)

Unconstrained K-MOD REV implant (UC, DC inserts):

Fixed inserts from 10 to 16mm. From K-MOD system:

- K-MOD TIBIA UC TRIALS (300110007)
- K-MOD TIBIA DC TRIALS (300110005)

Femoral implant too. Add the following trays:

- FEMORAL INSTRUMENTS (300113904)
- FEMORAL TRIAL BOXLESS (300113905)

Higher fixed insert from 19 to 25mm. Add:

- HIGH INSERTS TRIALS UC DC (300113906)

K-MOD REV CCK unconstrained (PS inserts) or constrained (CCK inserts) implant:

PS fixed inserts from 10 to 16mm from the K-MOD system:

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SAW BLADES

Resection slots do accept saw blades up to 1,47mm height.

The following oscillating saw blades are available:

- 120700050 K-MOD 1.47mm Oscillating saw Blade Stryker
- 120700051 K-MOD 1.47mm Oscillating saw Blade Linvatec MicroAire
- 120700052 K-MOD 1.47mm Oscillating saw blade Synthes Sodem

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
120700052
K-MOD 1,47mm Lama Oscillante attacco Synthes Sodem
K-MOD 1,47mm Oscillating saw blade Synthes Sodem
K-MOD 1,47 mm Hoja oscilante Synthes Sodem

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A narrow (approximately 13mm) oscillating blade or reciprocating saw blade is required for CCK femoral box preparation.

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PIN

The instrument set includes different types of pin. Unless specified otherwise, preferably use the Pin Dia 3,2 X 80mm (300110144) or Pin Dia 3,2mm L.120mm (300110278) depending on bone and guide positions.

Preferably insert the pin with the motor instead of the hammer, using the appropriate Jacobs Adaptor (300110377).

Always proceed by inserting the pin into the hole of the cutting guide by hand, even until contact with the bone, and then proceed to penetration into the bone, taking care to keep in line.

We suggest having around 6/8 pins available.

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Remove failed components and try to keep as much bone stock as possible.

Remove the PE insert and femoral and tibial components.


After cleaning the articulation, pre-check the existing flexion and extension gaps by using the appropriate Spacer Block REV FemPE 10/25mm (300113420/429), which is inserted on the EXT side to simulate the height of the femur and insert after having connected the tibial side Spacer Tibial Component REV (300113430) to simulate the tibial component.

If necessary, add the Augmentation Spacer Block REV 5mm (300113355). Remember that for each compartment a maximum of two (5+5=10mm) for the femoral side and three (5+5+5=15mm) for the distal femoral and tibial side may be added.

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
Following pre-op planning indications, progressively ream the IM canal up to the selected length and diameter using the sequential reamers Diaphysis Reamer Dia 10/20mm (300113150/160) connected to the T Handle (300113170). The last Diaphysis Reamer should just reach the cortical bone, avoiding press-fit of the stem into the canal.


Note: the length mark refers, for each stem length (60, 100, 140mm), to the maximum reachable depth including the offset module, whose length is equal to 20mm. The mark is also calculated on the femoral connection which is approximately 10mm longer than the tibial one.

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The last Diaphysis Reamer should be stable for the next steps. If necessary, fill the cavity of the failed tibial component with one Sleeve Small, Medium, Large (300113408/410) chosen depending by the cavity to be filled.

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Set the tibial cut as follows:

- Posterior tibial slope: 0°;
- Rotation: reference to third medial of the tibial tuberosity/center of the ankle; lock the IM Tibial Guide on the Reamer;
- Level of resection: use the Stylus REV 0mm/9mm (300113401) considering a refresh level (0mm) for revision and 9mm from the highest point of the condyle for primary surgery.

Fix the Cutting Guide REV with two pins in the 0 line.

Remove the IM Tibial Guide by first fully unscrewing the upper knob and then only slightly unscrewing the one on the side.

Unscrew the knob on the Diaphysis Reamer and remove it from the canal.

Stabilize the Cutting Guide REV with one or two pins in the oblique holes.

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
It's possible to check proper alignment by inserting the Tibial Baseplate Handle (300110196) with the Alignment Rod (300110161) in the Cutting Guide REV.

Perform tibial resection. If necessary, the level of the resection may be different in the medial and lateral parts (slots: 5, 10mm). If necessary, shift the Cutting Guide REV 2mm proximally or distally.

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

Choose the appropriate tibial size with the Baseplate REV 1/6 (300113201/206).

Based on chosen size and performed cuts, add one (5mm) or two (5+5=10mm) augmentations Tibial Augmentation 1/6, ML/LR - MR/LL 5mm (300113231/246) at the Baseplate REV, medial and/or lateral compartment. Evaluate a 15mm augmentation, medial and/or lateral compartment: Tibial Augmentation 1/6, ML/LR - MR/LL 15mm (300113251/266).

Note: 15mm augmentation borders develop a conical shape (5°) instead of cylindrical; as a consequence, carefully evaluate the match of the borders on the tibial resected surface and the need for a different size.

Place the last Diaphysis Reamer in the tibial IM canal; if necessary, fill the cavity with a Sleeve. Take care this is completely in place into the proximal IM canal; the proximal end should not affect the correct positioning of the Baseplate REV.

Evaluate if an off-set is necessary (3 or 6mm) or not (0), by using the appropriate Offset Tibial Bushing 0, 3, 6mm (300113220, 3, 6).

Assemble the Bushing with the Baseplate REV and slide on the Reamer down to the tibial cut. Rotational values (from 1 to 12) must be visible in the upper side. Rotate the assembly up to determine the correct tibial component position. In case of off-set, take note of the relative value (1 to 12).

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Fix the Baseplate REV with two pins in the anterior holes if possible, or in the ones on the upper side of the Baseplate. It's also possible to use the Headed Pin For Baseplate (300110197).
Remove the Bushing and Reamer.

Insert the Reaming And Broaching Keel Guide (300110199) by knocking the three spikes into the proximal tibia through the holes in the Baseplate REV.
Prepare the conical slot in the bone using the Conical Tibial Reamer (300110200) until it stops against the Keel Guide.

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Place the Stem Reamer Adaptor Tibia (300113230) into the Reaming Guide. Prepare the slot for the stem connection in the bone by using the Stem/Offset Reamer (300113402).
In case of tibial off-set (3, 6mm), the Stem/Offset Reamer should be introduced till the end. If not (0mm), it's enough to stop at the reference groove (see pictures).

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Finally, prepare the tibial keel slot using the appropriate size Keel Broach Size 1-2, 3-4, 5-6 (300110254/256), depending by the tibial size, assembled with the Tibial Broach Handle (300110271) and tapped till the end.

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Assemble the chosen trial components Trial Tibial Tray REV 1/6 (300113211/216) and Trial Stem 10x60/20x140 (300113100/132).

In case of off-set, between them assemble also the chosen Offset 3, 6mm Trial Adaptor (300113403, 406). Take care to position it in the pre-determined rotational position; the rotational value has to be positioned in line with the anterior mark on the connection of the tibial trial component.

Lock the assembly with the Trial Screw REV (300113407) using the Hex Screwdriver For Screws 4,5-6,5mm (340085045).


If used, under the Tibial Tray place the same Tibial Augmentations evaluated during the tibial preparation.

Add, in medial and/or lateral compartment, one (5mm) or two (5+5=10mm) Tibial Augmentation 1/6, ML/LR - MR/LL 5mm (300113231/246) or one Tibial Augmentation 1/6, ML/LR - MR/LL 15mm (300113251/266).

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Insert the trial assembly into the tibial canal. If necessary, slightly impact by using the Tibial Tray Impactor (300110258).

The tibial trial assembly should stay in place as a reference for the following femoral steps.

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Following pre-op planning indications, check femoral size and AP positioning by using the Femoral Template 1/2, 3/4, 5/6 (300110182/184).

Evaluate:

- AP and ML dimensions;
- Need for augmentations (indications: 5/10mm);
- Need for off-set module (indications: 3/6).

CCK version:

In the preliminary measurement phase, the Joint Line Ruler (300113375) allows you to evaluate, indicatively, the joint line and, therefore, the possible need for distal wedges.

Hook the Ruler into the hole on the Femoral Template of the chosen size.

Position the center of the Ruler hole in correspondence with the medial epicondyle; evaluate the possible need for distal wedges (5, 10, 15mm) on the graduated scale.

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
Following pre-op planning indications, progressive ream the IM canal up to the selected length and diameter using the sequential reamers Diaphysis Reamer Dia 10/20mm (300113150/160) connected to the T Handle (300113170). The last Diaphysis Reamer should just reach the cortical bone, avoiding press-fit of the stem into the canal.

The last Diaphysis Reamer should be stable for the next steps. If necessary, fill the cavity of the failed femoral component with one Sleeve Small, Medium, Large (300113408/410) chosen depending by the cavity to be filled.

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Length mark: next distal resection to be performed

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
300113408/410
Riempitivo Canale Small, Medium, Large
Sleeve Small, Medium, Large
Relleno Canal Small, Medium, Large

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Slide the assembly on the last Diaphysis Reamer.

Determine the level of resection using the Stylus REV 0mm/9mm (300113401) considering a refresh level for revision (0mm) and 9mm from the more prominent distal condyle for primary.

Fix the Cutting Guide REV with two pins in the 0 line.

Remove the whole assembly by first unscrewing the knob of the Distal Cut Guide then a little bit the one on the Reamer.

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
Stabilize the Cutting Guide REV with one or two oblique pins.

It's possible to check proper alignment by inserting the Tibial Baseplate Handle (300110196) with the Alignment Rod (300110161) into the Cutting Guide REV.

Perform the distal femoral resection; if necessary, the level of resection may be different medial and lateral (slot: 5, 10mm); if necessary, shift 2mm the Guide proximal or distal direction.

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
Check gap and ligament balancing in extension by using the Spacer Block REV FemPE 10/25mm (300113420/429), EXT side.

If appropriate distal cuts have been performed or if it's necessary to distalize the joint line, add one (5mm) or two (5+5=10mm) Augmentation Spacer Block 5mm (300113355) at the Spacer Block REV, distal femoral side, medial and/or lateral.


Caution: the check should be performed with the tibial trial already IN PLACE so the Spacer Tibial Component REV (300113430) should NOT be added to the Spacer Block REV.

If required, it's possible to use the Alignment Rod (300110161) in order to check mechanical alignment in proximal and distal directions.

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Based on the pre-determined femoral sizing, choose the appropriate Fem 4in1 Cut Block REV 1/6 (300113311/316).

Depending on the distal resection cuts and the extension gap check, eventually add on the back side of the 4in1 Cut Block, the appropriate Augmentation 5, 10, 15mm x 4in1 Cut Block (300113317, 318), 5 or 10 or 15mm, in the medial and/or lateral side.

Connect the prepared 4in1 Cut Block to the Femur Valgus Adapter 5° (300113302). Take care the correct leg indication (Left/Right) should be visible in the upper side.

Depending on pre-op plan and off-set evaluations during femoral sizing, connect also the appropriate Femoral Offset Bushing 0, 3, 6mm (300113340, 343, 346).

Note: the Offset Bushing may be changed during the next steps.

Place the last Diaphysis Reamer in the femoral IM canal; if necessary, fill the cavity with a Sleeve. Slide the assembly on the Reamer up to the distal femoral cut. Take care the Sleeve is completely in place into the distal IM canal; the distal end should not affect the correct positioning of the 4in1 Cut Block.

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In determining the 4in1 Cut Block final position, considered the following references:

- ANTERIOR REFERENCE

Set the anterior refresh or primary cut level (highest lateral side on the anterior femoral cortex) using the Stylus REV 0mm/9mm (300113401) in the anterior slot, 0mm reference.

- FLEXION GAP

Flexion gap should match the already created extension gap.

Pre-check flexion gap and set femoral extra-rotation by positioning the same Spacer Block REV FemPE 10/25mm (300113420/429) defined in extension, but now introducing the FLEX side between Trial Tibial Tray and the posterior side of the 4in1 Cut Block.

Determine the final position of the 4in1 Cut Block acting on:

- Femoral Offset Bushing 0, 3, 6mm: value and rotation, in order to modify AP and ML positions;
- Fem 4in1 Cut Block: rotation, in order to modify femoral component extra-rotation;
- Fem 4in1 Cut Block: size, in order to modify AP dimensions (posterior cut amount, flexion gap) and ML width.

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Once determined the desired position, lock the Femoral Offset Bushing in order to stabilize the assembly; fix the 4in1 Cut Block with two oblique pins, if possible with the Pin With Stop Diam 3,2mm L.40 (300110198).

Before to start resection, check the anterior cut with the Angel Wing (300110139).

In case of off-set, take note of the Bushing and the rotation value (from 1 to 12) as indicated on the Femur Valgus Adapter at the mark (Left, Right) corresponding to the operated leg.


Note: the control of the anterior resection is particularly indicated when using 6mm offset module and small femoral sizes, from 1 to 3.

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
In order to increase stability, it's possible to add two handles to the 4in1 Cut Block: the Universal Handle (300110279) and the Universal Handle/Pin Extractor (300110281).

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Unlock and remove the Femur Valgus Adapter and the Femoral Offset Bushing. Remove the Diaphysis Reamer from the IM canal.

Take care the two oblique pins are not in conflict with the Sleeve, not allowing the extraction. In case, slightly remove the pin, remove the Diaphysis Reamer and then fix again the pin.

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
Assemble to the 4in1 Cut Block: the Femur Valgus Adapter and the Stem Reamer Guide Femur (300113303). Prepare the slot for the stem connection in the bone by using the Stem/Offset Reamer (300113402). In order to increase stability, it's possible to add two handles also to the Stem Reamer Guide: the Universal Handle (300110279) and the Universal Handle/Pin Extractor (300110281). In case of femoral off-set (3, 6mm), the Stem/Offset Reamer should be introduced till the end. If not (0mm), it's enough to stop at the reference groove (see pictures). Take care to ream in line with the guide, avoiding tilting of the assembly.

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
If the CCK constriction is necessary, it will be necessary to proceed with the preparation of the central opening of the box using the 4in1 Cut Block for its correct centering with respect to the femoral component.

Remove the Stem Reamer Guide Femur. Press the button on the Guide Box 1/6 (300113372/4), chosen on the basis of the size and laterality (L/R) of the femoral component, and hook it to the Support Guide Box (300113371).
Connect the assembly to the 4in1 Cut Block and fix it with two 80mm Pins in the parallel holes, in the proximal and central rows.
Remove the oblique pins on the 4in1 Cut Block, press the button on the Guide Box and remove the Support and the 4in1 Cut Block.

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Slide the Box Guide until it contacts the anterior femoral bone, stabilize it with an oblique Pin 80mm in the most distal holes marked with laterality (L/R).




Carry out the center cut through the slots using an oscillating narrow or reciprocating saw.

Remove the pins and Box Guide and finalize the femoral box preparation

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Check gap and ligament balancing in flexion by using the same Spacer Block REV FemPE 10/25mm (300113420/429), determined in extension, EXT side.


If appropriate posterior cuts have been performed, add one (5mm) or two (5+5=10mm) Augmentation Spacer Block REV 5mm (300113355) at the Spacer Block REV, posterior femoral side, medial and/or lateral.

Caution: the check should be performed with the tibial trial already in place, so the Spacer Tibial Component REV (300113430) should NOT be added to the Spacer Block REV.

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Assemble the chosen trials components Femur Trial REV 1/6 Left/Right (300113321/336) or Femur Trial K-MOD REV CCK 1/6 Left/Right (300113381/396) and Trial Stem 10x60/20x140 (300113100/132). In case of off-set, between them assemble also the chosen Offset 3, 6mm Trial Adaptor (300113403, 406). Take care to position it in the pre-determined rotational position.

Lock the assembly using the Hex Screwdriver For Screws 4,5-6,5mm (340085045) and the screw:

- K-MOD REV: Trial Screw REV (300113407);
- K-MOD REV CCK: K-MOD REV CCK Screw Trial Implant (300113376).

If used, add on the Femur Trial, distal and/or posterior, medial and/or lateral, one (5mm) or two (5+5=10mm) Femoral Augmentation 5mm 1-2-3/4-5-6 (300113351, 352) chosen depending by femoral size.

In case of CCK femur implant it is possible to assemble up to three (5+5+5=15mm) distal femoral wedges in a medial and/or lateral position.

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300113403, 406
 Trial Adattatore Offset 3, 6mm
 Offset 3, 6mm Trial Adaptor
 Adaptor Prueba Offset 3, 6mm

300113321/336
 Trial Femorale REV 1/6 Sinistro/Destro
 Femur Trial REV 1/6 Left/Right
 Femur Prueba REV 1/6 Izquierda/Derecha

300113351, 352
 Spessore Femorale 5mm 1-2-3/4-5-6
 Femoral Augmentation 5mm Size 1-2-3/4-5-6
 Aumento Femoral 5mm 1-2-3/4-5-6

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
Insert the trial assembly into the femoral canal and impact using the Femoral Impactor (300110259).

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
CCK FEMUR sizes 1 and 2: If a 15mm distal thickness and a 10mm posterior thickness are used in the same compartment, the trial Femoral Augmentation should be positioned as follows:

- distal: flat part (marked 5mm) positioned towards the posterior;
- posterior: flat part (marked 5mm) positioned towards the distal.

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Select the appropriate fixed trial insert. Insert size is the same as femoral size. Insert height depends from the used Spacer Block REV.

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
FEMUR K-MOD REV CCK

Choose between inserts with Posterior Stability Trial Insert PS 1/6 x 10/16mm (300111301÷336) or constrained CCK Trial Insert 1/6 x 10/16mm (300113551/586).

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Remove the trials component using the Slide-Hammer (300110277).

Femoral trial assembly may be removed directly with the Slide-Hammer in the flat end.

Tibial trial assembly may be removed by screwing the Universal Handle (300110279) to the Trial Tibia and extracting with the Hook Extractor (300110380) connected to the Slide-Hammer.


Based on the trial reduction, open the final implants.

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Before to connect the stem, the tibial component and the off-set module (if required), accurately check the pre-assembled small headless locking screws in the posterior side of the tibial component keel and in the off-set module.

They do **NOT** have to be in conflict within the connections.

In case, un-screw them a little bit by using Allen Key L 2mm Grains (300113451).

In case of tibial off-set, position the appropriate module (3, 6mm) in the stem connection of the tibial component. Take care to position the module in the correct rotational position, using the anterior mark on the tibial keel as a reference.

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
Insert the stem into the appropriate diameter hole of the Stem Assembly Device (300113440); the connection'd be in the upper side. If the stem performs a diapason at the end, as a protection from the impaction, take care to slide into the diapason the Stem Diapason Protector (300113441).


Position the tibial component or the tibial component connected with the off-set module on the stem and impact just one time using a hammer and the Fixed Tibial Tray Impactor (300110258).

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

Safety lock the pre-assembled small headless locking screws in the posterior side of the tibial component keel and in the off-set module with the Allen Key L 2mm Grains (300113451).

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
In case of tibial augmentations, connect the appropriate wedge in the correct position medial and/or lateral after the removal of the relative plastic plugs. Lock the augmentation to the tibial plateau using the Allen Key L 3mm Wedge (300113450).

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Before to connect the stem, the stem connection and the off-set module (if required), accurately check the pre-assembled small headless locking screws in the lateral side of the femoral component

connection and in the off-set module.

They do NOT have to be in conflict within the connections.

In case, un-screw them a little bit by using Allen Key L 2mm Grains (300113451).

If no femoral off-set is required, directly proceed with step 2.

In case of femoral off-set, position the appropriate module (3, 6mm) in the stem connection of the femoral component. Take care to position the module in the correct rotational position, using the posterior mark on the femoral connection as a reference.

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
Insert the stem into the appropriate diameter hole of the Stem Assembly Device (300113440); the connection'd be in the upper side. If the stem performs a diapason at the end, as a protection from the impactation, take care to slide into the diapason the Stem Diapason Protector (300113441). Position the femoral component or the femoral component connected with the off-set module on the stem and impact just one time using a hammer and the Femoral Impactor (300110259).

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

Safety lock the pre-assembled small headless locking screws in the lateral side of the femoral component connection and in the off-set module with the Allen Key L 2mm Grains (300113451). If used, remove the Diapason Protector.

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In case of femoral augmentations, position the femoral component into the Femoral Augmentation Assembly Device (300113360).


Connect the appropriate wedges in the correct positions medial and/or lateral, distal and/or posterior. Lock the augmentations to the femoral component using the Allen Key L 3mm Wedge (300113450).

Suggestion: in order to facilitate the procedure, start with the distal wedges if used.

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
Accurately prepare the bone bed. Position bone cement in the back side of the femoral and tibial components. As a suggestion, the cement mantle should cover up to the stem connection or off-set module if used. Stems are intended for no cemented use.

Insert as much as possible the tibial component into the tibial canal and then impact with the Fixed Tibial Tray Impactor (300110258). Check for the correspondence between component and resections.

If bone screws are used, previously remove the plastic plug. Insert the cortical screw up the end taking cares:

- do not exceed 10° of inclination of the screw in the hole;
- the head of the screw does not overhang the upper surface of the tibial plateau in order do not conflict with the PE insert.

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Insert the definitive PE insert on the tibial component starting from the back side, light tilted; the groove on the back side should match the border of the tibial component. Lay down the anterior part too and lightly impact on the anterior border, oblique direction, up to the anterior connection.

Security lock the anti micromotion screw, included into the tibial component package, using the Hex Screwdriver For Screws 4,5-6,5mm (340085045).

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Flex the knee approx. more than 100°. Insert as much as possible the femoral component into the canal, then slightly extend the knee and impact with the Femoral Impactor (300110259). Check for the correspondence between component and resections.

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Accurately remove all the exceeding cement. Bring the leg in full extension up to when the cement is completely hard. Remove all the cement debris. Evaluate stability along ROM. Wash and clean the articulation. Close as surgeon's preferred procedure.

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In case of K-MOD REV CCK implant, in order to optimize the final femorotibial rotation, it is possible to proceed with the following alternative sequence.

TIBIA: RESECTION

Perform the resection without proceeding with the preparation and the trial implant.

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FEMUR: PREPARATION

Proceed with the complete femoral preparation. Note: Perform the external rotation setting in flexion using a baseplate instead of the trial tibia.

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SETTING THE TIBIAL ROTATION

Position the trial femoral component and engage, into its box, a CCK trial insert, alone or connected to a tibial baseplate.. Move the knee into flexion and extension and mark the correct rotation of the CCK insert to the tibia; this coupling will allow the correct final post/cam seating of the insert in the femur.

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TIBIA: PREPARATION

Proceed with the final preparation of the tibia using the mark made in the previous step as a rotational reference.

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
TRIAL REDUCTION FINAL IMPLANTS IMPLANTATION

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Optionally, with Trial Insert it's possible to perform a final check of PE insert parameters, as height and congruence type, even after final femoral and tibial components implantation.

FEMUR K-MOD REV

Choose between UltraCongruent inserts (UC): Trial UC Insert Size 1/6 Thkn 10/16mm (300110481/516) and Trial Insert UC 1/6 x 19mm (300113501/506) or medial pivot Dynamic Congruence (DC): Trial DC Insert Size 1/6 Thkn 10/16mm R/L (300110521/596) and Trial Insert DC 1/6 x 19mm R/L (300113511/526).

If a higher thickness is required (22, 25mm), assemble the Trial Insert 19mm to the Trial Spacer 3/6mm 1-2, 3-4, 5-6 (300113530/540) chosen size depending (1-2, 3-4, 5-6) and thickness to be added: +3 (=22mm) or +6 (=25mm).


Trial Spacer matches two sizes of Trial Insert (1-2; 3-4; 5-6).

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FEMUR K-MOD REV CCK

Choose between inserts with Posterior Stability Trial Insert PS 1/6 x 10/16mm (30011301÷336) or constrained CCK Trial Insert 1/6 x 10/16mm (300113551/586).

If a higher thickness is required (19, 22, 25mm), couple the Trial Insert 19mm to the Trial Spacer 3/6mm 1-2, 3-4, 5-6 (300113530/540) chosen on the basis of the size (1-2, 3-4, 5-6) and thickness to be added: +3 (=19mm) or +6 (=22mm) or +3+6 (=25mm).

The Trial Spacer fits pairs of Trial Tibial Insert sizes (1-2; 3-4; 5-6).

Connect the Trial Insert into the Trial Tibia and perform the trial reduction. Evaluate stability, joint line, ROM.

Take note of the final components of the prosthetic implant.

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

Temporary lock the Trial Insert on the definitive tibial component using the PURPLE Screw M5 Fixed Tibia (300110439).

In case of use of the Trial Spacer 3/6mm it will be necessary to fix the trial insert thus assembled using the M5 Trial Screw Inserts 22/25mm (300113542), YELLOW colour; to do this, use the hexagonal screwdriver for screw 4.5-6.5mm (340085045).

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After extension and flexion stability checks and once chosen the correct insert, remove the trial screw and insert.

Trial screws are colored in order to be different from the final one, included into the tibial component packaging.


Trial Screw M5 Fixed Tibia and Screw M5 insert 22/25mm do NOT have to be used for the final PE insert.

Proceed implanting the definitive PE insert as described in the previous chapter FINAL IMPLANTS IMPLANTATION, steps 2 and 4.

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
Remove any marginal osteophytes from the patella. If the thickness of the patella is to be measured, use the millimetric scale on the Patella Impaction Clamp (300110210).

Position the Stylus of the Patella Resection Clamp (300110205) at a resection level of around 8-10mm. Grip the patella in such a way that the stylus is in contact with the apex of the patella.

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
Evaluate the resection using the Angel Wing (300110139) in the lateral and medial slots and if necessary correct it by modifying the resection height of the stylus. Make particularly sure that the thickness of the residual patella is adequate and that the patellar ligament is secured.

Lock the Patella Resection Clamp, rotate the stylus and perform resection of the patella through the slots.

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
Measure the patella diameter with the Sizing/Drilling Guide Patella (300110185) and choose the appropriate size. The diameter of the patella component is independent of the size of the femoral component.

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
Drill three holes for the pegs of the patellar component using the Patella Peg Drill (300110189), drilling up to the stop. Alternatively, the Patella Impaction Clamp (300110210) may be used.

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
Position the Patella Trial Diameter 28, 32, 36 (300110186, 87, 88) of the chosen size on the patella and press lightly into the three holes. Check the patellar implant in the trial reduction.

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Carefully choose the diameter of the definitive implant and then proceed to apply the cement. Align the three fixation pegs with the three holes drilled in the patella and press lightly to insert the component.

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Complete the positioning with the Patella Impaction Clamp (300110210), whose screw adjustment makes it possible to keep the patellar implant compressed until the bone cement has set definitively.

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Make sure there is perfect adherence between component and resection and remove any excess cement.

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120701701/701
Patella Taglia 28, 32, 36mm
Patella Size 28, 32, 36mm
Patella Tam. 28, 32, 36mm

[OPZIONE: COMPONENTE ROTULEA INLAY]
[OPTION: INLAY PATELLAR COMPONENT]
[OPCIÓN: COMPONENTE PATELAR INLAY]

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

Select the diameter of the guide so that it is significantly within the patellar perimeter and center it on the central rim of the articular surface.

Coloque el K-MOD Módulo Plano Anterior Rótula (300110924) y el K-MOD Guía Rótula InLay D 20mm/26mm (300110919/21) en la K-MOD Pinza Rótula InLay Revestim (300110928).

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
Based on the chosen diameter, select the InLay Patella Reamer D 20mm/26mm (300110904/06) and hook it to the Universal Patella Reamer Rod (300110901).

Insert it into the Guide until it comes into contact with the patellar apex. Place the Reamer Calibration Thickness (300110903) on the upper edge of the Guide, selecting an appropriate milling thickness. Place the Patella Ring Stop (300110902) on the Shims and lock it in place.

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Remove the Cutter Calibration Shims and proceed with milling until the Stop Ring comes into contact with the upper edge of the guide.

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According to the chosen diameter, place the K-MOD InLay Patella Trial D 20mm/26mm (300110925/27) and perform the trial reduction.
Remove the Patella InLay Trial.

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

Select the definitive component based on the chosen diameter, proceed with the appropriate cementation, position it and fix it by tightening the forceps until the end of the cementation phase.

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
Fixed insert: we do not recommend the implant of a tibia of more than two sizes lower of the femur.

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Locking grain and M5 screw
for insert included

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[WEDGE TIBIALE 15mm ML/LR - MR/LL]
[TIBIAL WEDGE 15mm ML/LR - MR/LL]
[WEDGE TIBIAL 15mm ML/LR - MR/LL]

Ref.	DESCRIZIONE	DESCRIPTION	DESCRIPCIÓN
120702631	K-MOD REV Wedge Tib 15mm Tg.1 ML/LR	K-MOD REV Tib ML/LR Wedge Sz1 Th15	K-MOD REV Aumento Tib 15mm ML/LR Tm1
120702632	K-MOD REV Wedge Tib 15mm Tg.2 ML/LR	K-MOD REV Tib ML/LR Wedge Sz2 Th15	K-MOD REV Aumento Tib 15mm ML/LR Tm2
120702633	K-MOD REV Wedge Tib 15mm Tg.3 ML/LR	K-MOD REV Tib ML/LR Wedge Sz3 Th15	K-MOD REV Aumento Tib 15mm ML/LR Tm3
120702634	K-MOD REV Wedge Tib 15mm Tg.4 ML/LR	K-MOD REV Tib ML/LR Wedge Sz4 Th15	K-MOD REV Aumento Tib 15mm ML/LR Tm4
120702635	K-MOD REV Wedge Tib 15mm Tg.5 ML/LR	K-MOD REV Tib ML/LR Wedge Sz5 Th15	K-MOD REV Aumento Tib 15mm ML/LR Tm5
120702636	K-MOD REV Wedge Tib 15mm Tg.6 ML/LR	K-MOD REV Tib ML/LR Wedge Sz6 Th15	K-MOD REV Aumento Tib 15mm ML/LR Tm6
120702641	K-MOD REV Wedge Tib 15mm Tg.1 MR/LL	K-MOD REV Tib MR/LL Wedge Sz1 Th15	K-MOD REV Aumento Tib 15mm MR/LL Tm1
120702642	K-MOD REV Wedge Tib 15mm Tg.2 MR/LL	K-MOD REV Tib MR/LL Wedge Sz2 Th15	K-MOD REV Aumento Tib 15mm MR/LL Tm2
120702643	K-MOD REV Wedge Tib 15mm Tg.3 MR/LL	K-MOD REV Tib MR/LL Wedge Sz3 Th15	K-MOD REV Aumento Tib 15mm MR/LL Tm3
120702644	K-MOD REV Wedge Tib 15mm Tg.4 MR/LL	K-MOD REV Tib MR/LL Wedge Sz4 Th15	K-MOD REV Aumento Tib 15mm MR/LL Tm4
120702645	K-MOD REV Wedge Tib 15mm Tg.5 MR/LL	K-MOD REV Tib MR/LL Wedge Sz5 Th15	K-MOD REV Aumento Tib 15mm MR/LL Tm5
120702646	K-MOD REV Wedge Tib 15mm Tg.6 MR/LL	K-MOD REV Tib MR/LL Wedge Sz6 Th15	K-MOD REV Aumento Tib 15mm MR/LL Tm6

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




















^{EN}
Use the QR-Code to view the instruction for use

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The informations in this document refers to the basic system configuration and sizes.

To check the complete range, refer to the relative technical data sheet.

Some codes / products on request, developed for specific markets, may not be in immediate availability.

Check availability with your local GRUPPO BIOIMPIANTI distributor.

Enjoy Mobility

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Via Liguria 28, 20068 Peschiera Borromeo (Milan) Italy

Tel. +39 02 51650371 - Fax +39 02 51650393

info@bioimpianti.it

bioimpianti.it