Enjoy Mobility

KINEMATIC ALIGNMENT

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Following the introduction of the new K-MOD TOUCH 3REF instrumentation, which enables the use of an anterior, posterior, or hybrid femoral reference, the K-MOD system has undergone further development and now includes the new additional K-MOD KA series of instruments, which makes it possible to use the kinematic alignment technique for a total knee replacement.

This technique aims to recreate the articular kinematics of the patient's native knee; the technique can therefore be considered specific to the individual patient. Kinematic alignment is an option available to orthopaedic surgeons to minimise the risk of pain arising from a total knee replacement. Indeed, many scientific studies show that a significant percentage of patients who undergo mechanical alignment are still dissatisfied due to pain, instability or stiffness in the joint.

Kinematic alignment aims to work according to the functional phenotype of the individual patient, aligning the knee prosthesis according to their specific anatomy.

Kinematic alignment is therefore based on the following principles:

1. Restoration of the native alignment of the patient, knee and leg
2. Restoration of the native femoral-tibial joint lines of the knee
3. Restoration of the native laxity of the knee

The kinematic alignment technique therefore minimises the need for capsular, collateral, or posterior cruciate ligament release, and is essentially based on bone resection. The importance of this is highlighted by the caliper, which is used, for example, to measure the distal and posterior femoral resection in order to restore the bone using a prosthesis with an identical thickness.

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Usura cartilagine: 0mm

Spessore lama: 1mm

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The K-MOD TOUCH 3REF and K-MOD KA instruments include:

- Intramedullary femoral guide for distal resection with the possibility of separately compensating 1 or 2 mm of worn medial and lateral articular cartilage - External rotation guide for femur, with posterior reference
- Caliper for measuring distal and posterior femoral bone resection
- Cartier tibial resection guide for 3° varus and stylus for adjustable resection height of 0-10mm
- Modular spacer blocks of 10 to 16mm, to set the tibial resection based on the creation of a uniform rectangular gap in extension and/or flexion
- 2° varus/valgus tibial correction block
- 1 to 4° femoral extrarotation correction block

For correct kinematic positioning of the femoral component, the distal and posterior femoral resections must be equal to the thickness of the prosthesis: 9mm. This assessment includes the possible compensation, through suitable thicknesses, of the wear of the articular cartilage (1 or 2mm), and the thickness of the oscillating blade used for the resections is assumed to be 1mm.

The thickness of the femoral component, whether distal or posterior, medial or lateral, is 9mm, and must be determined from the sum of the following: bone resection + cartilage wear (0, 1, 2mm) + thickness of blade (1mm).

The accuracy of the resection is checked by measurement using the caliper.

Example of medial distal condyle

Cartilage wear (compensated by distal thickness): 2mm

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On the patient's most recent X-rays, perform thorough pre-op planning on the AP and ML planes by using the X-Ray templates K-MOD: K-MOD X-Ray Template Fem AP/ML, Tib AP/ML, Fix/Mob, Stem (KMOD-XR-000÷008; KMODPS-XR-000÷003).

Evaluate (magnification: 15%): sizes of the femoral and tibial components, joint line, diameter and length of the optional tibial stem if needed.

Evaluation is indicative and must always be checked and modified during surgery if necessary.



Technologies are available for personalized 3D preoperative planning through online platforms.

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
The KA kinematic alignment technique involves use of the basic K-MOD TOUCH 3REF instrument set consisting of basic trays:

- TIBIAL INSTRUMENTS (300110003);
- 3REF FEMORAL INSTRUMENTS (300200015);
- FEMORAL TRIAL (300110002);

Trial inserts:

- K-MOD TIBIA CR TRIALS (300110007);
- K-MOD TIBIA UC TRIALS (300110007);
- K-MOD TIBIA DC TRIALS (300110005);

SAW BLADES

In the kinematic alignment technique, the precision of resection is essential to the success of bone preparations. For this purpose, it is essential to choose fairly rigid blades of appropriate thickness. Resection slots do accept saw blades up to 1,47mm height. No need for reciprocating blades.

The following oscillating saw blades are available:

- 120700050 K-MOD 1,47mm Oscillating saw Blade Stryker
- 120700051 K-MOD 1,47mm Oscillating saw Blade Linvatec MicroAire
- 120700052 K-MOD 1,47mm Oscillating saw blade Synthes Sodem

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
120700052
K-MOD 1,47mm Lama Oscillante attacco Synthes Sodem
K-MOD 1,47mm Oscillating saw blade Synthes Sodem
K-MOD 1,47 mm Hoja oscilante Synthes Sodem

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
Assess any wear on the medial and lateral distal condylar surfaces; the system can compensate for 1 or 2mm of missing articular cartilage, independently in the two condyles. In the event of wear, smooth out the affected condylar surface appropriately.

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
Use the Femoral IM Reamer 8.25 mm (300110145) to prepare the intramedullary hole in the femoral canal. The drill hole must be located approximately 1 cm above the femoral attachment of the posterior cruciate ligament and slightly medial to the midline of the femur. Removal of osteophytes from the margins of the intercondylar notch may aid identification of landmarks. Insert the SHORT IM Rod (300110376) into the femoral canal hole using the relative T-Handle For IM Rod (300110151). The IM Rod has longitudinal grooves that limit pressurization of the canal and the risk of emboli. Remove the T-Handle

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Assemble the K-MOD 3REF Femoral Sizer ExtraRot (300110702), the K-MOD 3REF Femoral Valgus Guide (300110701) with the word FRONT facing upwards and, in the event of wear, select the corresponding K-MOD KA Distal Spacer 7/8mm Cut, (300110730, 31) to be attached to the back of the Sizer at the site of the worn condyle. The Distal Thickness for Cuts can be used to compensate for wear, specifically:

8mm CUT - Wear: 1mm

7mm CUT - Wear: 2mm

In the example, the Distal Thickness only compensates for medial wear, estimated at 2mm.

Insert the distal femoral resection guide onto the IM Rod until the ExtRot Femur Sizer is in contact with at least one of the two femoral condyles.

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


Position the compass of the Femoral Valgus Guide so that the marking relating to the working leg, L (Left) or R (Right), is visible on the side facing upwards. Rotate it until both the medial and lateral sides of the Femur Sizer (or Thickness of Distal Cuts where applicable) are in contact with the respective condyles.

This will enable resection of the correct thickness of distal bone independently in the two distal condyles, which will be 9mm (if there is no wear) or 8mm (if there is 1mm of wear) or 7mm (if 2mm of wear), then accurately reconstruct the native distal femoral joint line, since the thickness of the femoral component is 9mm.

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
Finally, assemble the Distal Femur Cutting Block (300110159) so that the wording DISTAL FEMUR can be seen facing upwards, and secure it using two 80mm pins in the parallel rows of holes marked 0.

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
Release the Guides by pressing the button on the side of the Distal Femur Cutting Block and slide them off along the IM Rod. Remove the IM Rod from the femoral canal by the T Handle. The stability of the Distal Femur Cutting Block can be increased by inserting one or two Pins into the slanted holes of the Distal Femur Cutting Block. Perform the distal femoral resection.

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In the example, the distal thickness of the femoral component (9mm) reconstructs the native distal femoral joint line as:

Medial condyle: Wear (2mm) + Blade thickness (1mm) + Resected bone (6mm)

Lateral condyle: Wear (0mm) + Blade thickness (1mm) + Resected bone (8mm)

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



Note: the thickness of the blade is assumed to be 1mm.

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Assemble the K-MOD 3REF Femoral Sizer Extra-Rot (300110702) and the Fem ExtraRotation Guide (300110703). It is assumed that the posterior condyles are not subject to significant wear. Therefore in order to resect an amount of bone that is equal to the thickness of the posterior femoral component (9mm), select a femoral extra-rotation value of 0°. Loosen the clamping screw with the screwdriver (340085045); select this value by bringing the arrowhead towards it, then tighten the clamping screw.


With the leg at about 90° flexion, put the assembled components into position and make sure that the Femoral Sizer Extra-Rot is completely in contact with the distal femoral cutting plane and that the adjustable Fem ExtraRotation Guide is in contact with the posterior condyles. Set the Femoral Stylus by sliding it into its site; the SIZE CHECK indication must be in correspondence to the value of the chosen measure during preoperative planning.

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Insert two 3.2mm x 80mm Pins (300110144) into the upper holes of the 3REF Femur Sizer indicated with the letter P* (Posterior Reference).

To check the position of the anterior resection, insert the Angel Wing into the upper, lateral or medial slot of the 3REF Femur Sizer.

Place the tip of the Femoral Stylus (300110167) on the anterior cortex towards the highest lateral part in order to avoid the problem of anterior femoral notching.

Read the femoral size through the front slot of the Femur Sizer; if it does not correspond to the stylus setting, reset the stylus and and recheck the size.

If the indicator is midway between two sizes, the surgeon may choose the smaller or the bigger.

In the POSTERIOR reference TOUCH instrumentation procedure, the posterior resections are the same for each size of the femoral component and, in neutral rotation, equal to the posterior size of the femoral component (9mm).

By choosing the larger size, a larger portion of anterior bone is preserved, the risk of notching is reduced, but, there is the risk of excessive tension in the patellofemoral joint.

By choosing the smaller size, a smaller portion of anterior bone is resected; the possibility of excessive tension in the patellofemoral joint is reduced, but there is a risk of creating a femoral notching. In order to limit this risk, the anterior femoral resection develops a 4° inclination.

Rotate the stylus on the medial side and remove the assembly from the femur.

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

The 3REF Resection Block 4IN1 (300110711÷6) has an outside measurement in the ML plane equal to that of the relative femoral component and it also shows the marks relating to the measurement in ML also for the Skinny version (sizes 2-, 3-, 4-).

Position the 3REF Resection Block 4IN1 corresponding to the chosen measurement on the Pins, using the 0 holes on the bottom row, associated with the posterior reference (POST REF). Slide it until it fully adheres to the distal resection plane of the femur.

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Check the position of the anterior resection with the help of the Angel Wing. In order to avoid the problem of anterior femoral notching, the exit angle of the anterior resection should be equal to 4°.

The level of anterior resection can be varied by changing the size of the femoral component. For example, in case of danger of notching, it is possible to evaluate a higher measure that shifts the resection level towards the anterior. The difference in AP between two different femoral sizes is 4mm, while between size 5 and 6 it is 5mm.

Secure the Resection Block 3REF 4IN1 into its final position by inserting two 40mm short pins into the slanted holes.

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

The posterior thickness of the femoral component (9mm) reconstructs the native posterior femoral joint line as:

Medial condyle: Wear (0mm) + Blade thickness (1mm) + Resected bone (8mm)

Lateral condyle: Wear (0mm) + Blade thickness (1mm) + Resected bone (8mm).


Use the Caliper (300110206) to assess and confirm the exact extent of the medial and lateral condylar resections.

Note: the thickness of the blade is assumed to be 1 mm.

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

Perform the anterior resection; remove the parallel pins and proceed with the posterior and anterior oblique cuts. Remove the Pins and the 3REF Resection Block 4IN1

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Assemble the extramedullary tibial resection guide: the Spikeless Tibial Resection Proximal Guide (300110115), the appropriate Tibial Cutting Block RIGHT/LEFT (LEFT 300110741, RIGHT: 300110742), the Tibial Distal Guide (300110122) and the Ankle Clamp (300110133). Raise the ferrule of the Spikeless Tibial Resection Proximal Guide so that its slot is in line with that of Tibial Cutting Block RIGHT/LEFT. Secure the malleoli with the Ankle Clamp and adjust the length of the guide so that the Tibial Cutting Block RIGHT/LEFT approximately at the height of the desired cut. Align the tibial resection guide on the proximal tibia with the medial third of the tibial tuberosity; in the sagittal plane, the guide must be parallel to the anterior tibial cortex. Lock it in place by passing an 80mm temporary Pin through the slot of the Tibial Cutting Block RIGHT/LEFT; the temporary Pin will still allow any necessary correction of the posterior slope, the varus-valgus mechanical axis and the resection height. The posterior slope may be adjusted by using the Angel Wing (300110139) and operating the knob on the Spikeless Tibial Resection Proximal Guide and/or the AP trigger on the Tibial Distal Guide. In general, proximal tibial resection should follow the tibia's natural posterior slope (medial condyle). It is advisable not to exceed the posterior slope.

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The alignment of the frontal mechanical axis in the varus-valgus plane may be adjusted by unscrewing the knob on the Tibial Distal Guide and sliding the guide on the Ankle Clamp in the ML plane

For correct varus-valgus alignment, please refer to the following sections: TIBIA: SETTING WITH GAP IN EXTENSION and TIBIA: SETTING WITH GAP IN FLEXION.


Position the K-MOD KA Tibial Stylus Adjustable (300110743) in the slot on the Cutting Block. The Stylus can be adjusted for resection from 0 to 10mm. The surgeon can choose to use the lowest point of most worn side of the tibial plateau, or the highest point of the least worn side, as the reference for the cut.

Position the tip of the Stylus on the chosen bone reference, then lower the ferrule of the Spikeless Guide so that the cutting height can be adjusted. Turn the Stylus handgrip to determine the most suitable resection, then reposition the ferrule in contact with the Cutting Block.

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
It is possible to finalise the tibial resection setting, and check the extension gap, in advance. Assemble D/P Fem Modular Reference (300110735) with a K-MOD KA Modular Spacer Tibia L/R 10÷16 (300110750÷59) chosen according to the side involved (LEFT/RIGHT) and the width of the gap (10, 12, 14, 16mm).

With the Tibial Guide in place, place the Distal Femur Cutting Block I (300110159) in the chosen position for distal resection. Bring the knee into full extension and insert the D/P Fem Modular Reference paddles into the distal femoral resection slot with the K-MOD KA Modular Spacer Tibia L/R 10÷16 facing the tibia.

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
Adjust the resection height by turning the ferule micrometrically, and the tibial varo-valgus alignment by releasing the handgrip of the Tibial Distal Guide and sliding the guide on the Ankle Clamp in the ML plane until the lower surface of the Modular Spacer is fully in contact with the upper surface of the Tibial Resection Block.

This adjustment makes it possible to set a tibial resection parallel to the distal femoral resection and thereby create a rectangular extension gap whose width is equal to the thickness of the K-MOD KA Modular Spacer.

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


Secure the Tibial Cutting Block using two in the parallel rows of holes marked 0. The stability of the guide can be increased by inserting a third pin into the slanted holes on the medial side of the Tibial Cutting Block.

Remove the femoral assembly and bring the knee into flexion. Check that the lever of the Tibial Cutting Block is correctly closed and resect the proximal tibia through the appropriate slot, keeping the extramedullary tibial resection guide in place.

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It is possible to finalise the tibial resection setting, and check the flexion gap, in advance. Assemble the K-MOD KA D/P Fem Modular Reference (300110735) with a K-MOD KA Modular Spacer Tibia L/R 10÷16 (300110750÷59) chosen according to the side involved (LEFT/RIGHT) and the width of the gap (10, 12, 14, 16mm).



With the Tibial Guide in place, place the 3REF Resection Block 4IN1 (300110711÷6) in the chosen position for femoral preparation.

Bring the knee into 90 degree flexion and insert the paddles of the K-MOD KA D/P Fem Modular Reference into the slits for posterior femoral resection using the K-MOD KA Modular Spacer Tibia L/R 10÷16 facing the tibia.

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
Adjust the resection height by turning the ferrule micrometrically, and the varo-valgus tibial alignment by releasing the handgrip of the Tibial Distal and sliding the guide on the Ankle Clamp in the ML plane until the lower surface of the Modular Spacer Tibia L/R 10÷16 is fully in contact with the upper surface of the Tibial Cutting Block.

This adjustment makes it possible to set a tibial resection parallel to the posterior femoral resection and thereby create a rectangular flexion gap whose width is equal to the thickness of the Modular Spacer.

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
Secure the Tibial Cutting Block using two Pin 80mm in the parallel rows of holes marked 0. The stability of the guide can be increased by inserting a third pin into the slanted holes on the medial side of the Tibial Cutting Block.

Remove the femoral assembly and finish the femoral resections if not already performed. Check that the lever of the Tibial Cutting Block is correctly closed and resect the proximal tibia through the appropriate slot, keeping the extra-medullary tibial resection guide in place.

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

If necessary, the frontal alignment of the tibial resection, whether set up or already performed, can be adjusted using the Tibial Cutting Block. The Angular 2° Cutting Block (300110166) actually enables an adjustment of 2 degrees, in varus or valgus. In addition to this, it enables this adjustment by increasing the thickness of the resection by 2 or 4mm compared with the original (black line).

Insert the Angular 2° Cutting Block onto the same parallel pins used in the original tibial resection, using the same line of holes used for the Resection Block (0, +2, +4). If an angular correction and an increase in the depth of resection from the 0 line are required, use the +2 or +4 lines.

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
Depending on the side used, a correction of 2° in valgus (green line) or 2° in varus (red line) can be made.


The stability of the Angular 2° Cutting Block can be increased by inserting a third Pin into the slanted holes on the sides. Resect the proximal tibia through the slot.

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If necessary, the extrarotation of the femoral component can be adjusted after inserting the two pins into the K-MOD 3REF Femoral Sizer ExtraRot. The Fem ExtraRot Change Block (300110736) makes it possible to change the alignment from 1 to 4 degrees, in intra or extra rotation.

Remove the 3REF Resection Block 4IN1. Select the side of the Fem ExtraRot Change Block to be used according to the direction of the central arrow, which shows the direction of rotation with respect to the original alignment (black line).


Insert the Fem ExtraRot Change Block onto the Pins through the holes indicated by the arrows; to be specific, use the upper holes to change the rotation by 1 or 2 degrees or the lower holes to change the rotation by 3 or 4 degrees.

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
Secure the Fem ExtraRot Change Block using two Pins in the oblique holes at its ends. Remove the two parallel Pins and insert them into the two holes showing the value of the rotation with respect to the original (1/1, 2/2, 3/3, 4/4). Remove the two Pins in the oblique holes.

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

Remove the Fem ExtraRot Change Block. The two parallel Pins are now oriented according to the newly chosen femoral rotation (red line) as opposed to the original one (black line). Insert the 3REF Resection Block 4IN1 into the parallel Pins, still using the row of holes previously used (0, +/-1, +/-2); secure it through two Pins into the slanted holes at the ends and proceed with the AP and oblique resections.

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KA Operating technique video

EN

Use the QR-Code to display the video that simulates the operating technique

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IFU

EN

Use the QR-Code to view the instruction for use

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The informations in this document refers to the basic system configuration and sizes.

To check the complete range, refer to the relative technical data sheet.

Some codes / products on request, developed for specific markets, may not be in immediate availability.

Check availability with your local GRUPPO BIOIMPIANTI distributor.

Enjoy Mobility

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Via Liguria 28, 20068 Peschiera Borromeo (Milan) Italy

Tel. +39 02 51650371 - Fax +39 02 51650393

info@bioimpianti.it

infofrance@bioimpianti.it

bioimpianti.it