

**FOLLOW UP 1 MONTH AFTER THE FIRST CASE AND IMPLANT OF K-MOD  
FOLLOWING A UNI REVISION, WITH THE NEW INSTRUMENT SET  
K-MOD UNI REV**


Surgeon: Dr. Carmine Naccari

Hospital: Istituto di Cura Città di Pavia, Pavia








**EXTRACTOR AND CHISELS**  
to extract the failed  
unicompartmental prosthesis



**CALIPER**  
to measure the thickness of removed component and  
choose the distal and posterior wedges





**DISTAL AND POSTERIOR WEDGES**  
Fixation to K-MOD 3REF Femoral Sizer  
Thicknesses of distal and posterior wedges: 4, 5, 6, 7,  
8, 9mm